

Date: Wednesday, April 13, 2011
To: Lumenera Scientific Partners
From: Darren Bessette, Product Manager



Re: **INFINITY1-2 Product Modification Notice**
New Image Sensor Due To Manufacturer Obsolescence
Product News #11-Sci-02

Due to a manufacturer's part obsolescence, Lumenera will now begin phasing in a new version of its INFINITY1-2 2.0 megapixel color USB 2.0 camera that will include an updated sensor.

Although some camera specifications have been slightly modified, overall you should expect similar, high-quality performance to the legacy INFINITY1-2. This product modification will not affect pricing.

Please note that as a result of this modification:

- By default the camera is set to run at 10 frames per second (fps), with the option of increasing 15 fps
- New camera drivers are required - the camera software included with purchase contains both the updated drivers as well as the older versions

Any open orders will be transitioned to the new camera model once the current stock is depleted. You will be able to distinguish the new version by updated part numbers on the camera label and the change in external power from 6V to 5V (please note that external power is not typically required as it is powered by USB).

New ordering information as of April 13, 2011:

Product names and ordering part numbers have been modified as follows:

Camera Description	Previous Part Number	New Part Number
2.0 megapixel enclosed color camera	INFINITY1-2C	INFINITY1-2CB

Customers will continue to receive existing levels of technical support for these cameras as per our standard practice. Standard warranties apply. Lumenera reserves the right to substitute new parts as required on products returned to the factory for repairs. This policy may be subject to other customer-specific service level agreements.

We apologize for any inconvenience this may cause your organization. Our support and sales staff are available to assist you in any way possible. To address any questions regarding order, pricing and availability please contact your Business Development Representative or Isabel Martins, Sales Coordinator, isabel.martins@lumenera.com / +1-613-736-4077 x 180. For Technical Support please contact us at support@lumenera.com.

We thank you for your cooperation in this product transition.