

Date: Tuesday, July 17, 2012
To: Lumenera Scientific Partners
From: Darren Bessette, Product Manager



Re: Product Modification Notice
INFINITY USB 2.0 Enclosure Modifications
Product News #12-SCI-04

In order to improve functionality, over the coming months Lumenera will be phasing in two different modifications to INFINITY USB 2.0 camera enclosures:

Tripod Mount Modification:

An indentation has now been carved next to the tripod mount, providing a flat, stable surface for the camera to sit completely perpendicular when mounted on a tripod. Included are small spacers that can be used around the 1/4-20" screw, allowing a more secure surface for mounting. Please refer to the Application Note entitled "Alternative Mounting Options for INFINITY USB 2.0 Cameras" for additional information.



A package of two spacers can be ordered using the Lumenera part number La051.

Lens Mount Alterations:

Changes to the lens mount are also being phased in. The locking mechanism has been modified in order to provide a more secure camera-microscope connection and simplified detachment.



Both of these changes will apply to all INFINITY1, 2, X and HD cameras, as well as the INFINITY3-1U. The INFINITY3-1 and INFINITY4-11 will not be affected.

Please note: When mounting the camera to a microscope port, the spacers are not required.

To address any questions or concerns please contact your Business Development Representative or Isabel Martins, Senior Sales Coordinator, at isabel.martins@lumenera.com / +1-613-736-4077 x 187. For support please contact support@lumenera.com.